LIFEVIT

BÁSCULA INTELIGENTE
ESCALA INTELIGENTE
SMART SCALE
PÈSE-PERSONNE INTELLIGENT

KRYOS

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ENGLISH

IMPORTANT WARNING

**Read the instructions carefully
before using the product. Keep them
in a safe place for future reference**

SAFETY MEASURES

- This device is intended solely for domestic use.
- Remove the batteries from the product if it is not being used for a long time.
- Keep this appliance out of the reach of children and people with disabilities.
- Do not move the device while in use.
- Do not exceed the weighing capacity of the device.
- Do not use the device with wet feet.
- Do not store the device in upright position.
- Do not use the device with cracked or broken glass.
- Improper use of the product may involve danger and revoke the warranty and liability

Lifevit accepts no liability for any damage to people, animals or objects as a result of not observing these warnings.

GENERAL INFORMATION:

We would like to thank you for your trust in Lifevit

First check that the appliance is complete and undamaged. Contact the store where you purchased the product in the event of any doubt.

With your purchase you have acquired:

- 1 Lifevit scale, battery operated (4xAAA), maximum weight 180 kg, resolution 100 g. - batteries not included

INSTRUCTION FOR USE

The safety instructions marked at the beginning of this manual must be followed at all times when using this product.

Follow the instructions below in order to switch on the device:

1. Remove the product from its original packaging taking special care not to damage any of its parts. Do not use excessive force to remove the product from the box.
2. Remove the cover of the battery on the back of the scale.
3. Place the charged battery in the right place under the cover following the indicated polarity, that is, placing the positive "+" facing outwards.
4. Put the cover back on.
5. Place the device on a hard, rigid, flat, stable horizontal surface free of vibration.

NORMAL WEIGHT MODE

1. Gently place your feet on the surface of the scale to start the device and wait for the symbol "0.0" to appear on the screen.
2. Place yourself in the centre of the platform and avoid sudden movements, keep still and without touching other objects, until the scale shows your weight. The user must wait while the weight is being calculated. The flashing final value indicates that the weighing process is completed.
3. To save energy, the device goes into sleep mode after several seconds of inactivity.
4. To return to normal operation, simply place your foot lightly on the platform and wait for the symbol "0.0" to appear on the screen.
5. Replace the unit in a safe and stable place.

WEIGHT MODE WITH COMPOSITION INFORMATION

1. This product enables you to obtain information regarding the composition weight, indicating the values for body mass index, bone, muscle, water, etc. To access this mode, use the Lifevit app (see section Connectivity with Lifevit app)
2. Users enter their personal details via the Lifevit app, indicating their age, height and gender.
3. Once the user has been saved in the Lifevit application, there appears a link to the scale in the main window of the application. Press on the Bluetooth connection symbol in the upper right corner of the window while switching on the scale. The symbol turns blue when the devices are connected and ready to measure.
4. Gently place your bare feet on the surface of the scale to start the device. You should ensure that your feet are in contact with the four squares on the product surface.
5. Place yourself in the centre of the platform and avoid sudden movements, keep still and without touching other objects, until the scale shows your weight. The user must wait while the weight is being calculated.
6. It will then proceed to calculate the composition of your weight, which takes about

5 seconds. The values of BMI, fat, water, muscle, bone, Kcal will be displayed on screen after this calculation. Calories are shown on screen referring to the basal metabolic rate, i.e. the number of daily calories needed to maintain the basic functions of your body.

7. You may get off the device only after that calculation of the composition is performed.

8. The data will be displayed twice and then the scale will turn off.

Connection with Lifevit app:

This product can be connected with Lifevit app, downloadable on your Android or IOS device.

To connect the product to the app, go to the home screen of the app, proceed first to select the corresponding user in the mobile application. Once selected, you will be able to link the scale with the mobile device by pressing the Bluetooth symbol in the upper right corner on the main screen of app at the same time as the scale is switched on.

When the scale is connected with the mobile device the Bluetooth symbol will turn blue. Proceed to perform the weight measurement and manage your data easily in the Lifevit app.

Contraindications:

This product is for domestic use only and not for professional purposes.

Children must be supervised by an adult to prevent them using the device as a toy.

This device is not designed to be used by persons (including children) with reduced physical, sensory or psychological capacities, or who lack the necessary experience and knowledge, unless they are supervised by a person responsible for their safety. Contact your doctor if you are pregnant, or have a pacemaker, artificial limbs, electronic implants.

Repair:

Only specialized personnel should perform repairs to this electrical appliance. Any incorrectly performed repair can lead to considerable hazards for the user. In the event of repair, please always contact our aftersales service. Failure to observe these instructions will lead to the total loss of the product warranty.

TROUBLESHOOTING GUIDE: Possible operating anomalies

INCIDENT	CAUSE	SOLUTION
The device does not work	Device not switched on or badly positioned battery	Put on the device, and wait until your weight is shown or press on it with your feet until "0.0" is shown on the display. Check that the battery is installed with positive polarity "+" facing outwards.
	The device is damaged	Contact the Lifevit technical service.
Error message "Err"	Overload	It means that the maximum weight is exceeded.

PRODUCT CLEANING

This product can only be cleaned by hand, using a soft damp cloth and mild soap. The use of abrasive cleaning products or brushes, solvents, benzine, thinner or alcohol is completely prohibited. Never immerse the device or any of its parts in water or never moisten it excessively in order to clean it. Do not re-use the device until it is completely dry.

PRODUCT DISPOSAL

Do not dispose of the product in the litter at the end of its shelf-life. Contact the organic collection centre nearest to your home in accordance with the Waste Electrical and Electronic Equipment (WEEE) Directive 2012/19/EU to contribute to the protection of the environment.

Request detailed information regarding the current methods for disposing of packaging and used devices from your distributor.

We hope that you will be satisfied with this product.

SYMBOLS

This product complies with European Directive 2012/19/EU on electrical and electronic devices, known as WEEE (Waste Electrical and Electronic Equipment), that provides the legal framework applicable in the European Union for the disposal and reuse of waste electronic and electrical devices.

Do not dispose of this product in the litter. Go to the electrical and electronic waste collection centre closest to your home.

The product may contain batteries. Remove them before disposing of the product and dispose in special containers approved for this purpose.

This symbol means you should be careful and that the product surface is slippery when wet.

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WARRANTY REPORT

Lifevit. guarantees compliance of this product for the purpose for which it is intended for a period of two years for the entire Spanish territory. In the case of breakdown during the term of this warranty, you are entitled to repair or to the replacement of the product at no charge if the former is unavailable, unless one of these options prove impossible or disproportionate. In this case, you can then choose a reduction in price or cancellation of the sale, which must be dealt with directly with the sales vendor. This also covers replacement of spare parts provided that the product has been used according to the recommendations specified in this manual for both cases, and has not been tampered with by any third party that is unaccredited by Lifevit. This warranty does not affect your rights as a consumer in accordance with the provisions of Law 23/2003 dated 10 July 2003 on warranties on the sale of consumer goods in the Spanish territory.

USE OF WARRANTY

Consumers should contact a technical service centre authorised by Lifevit. for repair of the product as any tampering thereof by persons not authorised by Lifevit., as well as misuse or inappropriate use, will lead to the cancellation of this warranty. The warranty must be fully completed and delivered along with the receipt or delivery docket for the effective exercise of rights under this warranty.

Claims must take place within three years from delivery of the product according to the provisions of the warranty, and the consumer must inform the sales vendor of the lack of conformity, provided that such claim is made within the statutory period of warranty. This warranty should be retained by the user as well as the invoice, receipt or the delivery docket to facilitate the exercise of these rights. For technical service and after-sales care outside the Spanish territory, please submit your query to the point of sale where you purchased the item or by contacting Lifevit.

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